

# GoCociety®



GoLive  
DATABASE RECORDS v1.2

# GoLive database records

This document identifies all information records gathered and stored in the GoLive environment, regarding the GoLiveClip, GoLivePhone and GoLiveAssist components. To allow a better understanding the information is segmented in modules:

## Emergency Contacts

- Email (main, 1st optional and 2nd optional)
- Name (main, 1st optional and 2nd optional)
- Phone (main, 1st optional and 2nd optional)
- Use call center
- Use cascading alarms

## Call Centers

- Logo
- Name
- Phone number

## Device Info

- Android Version
- Android ID
- GoLivePhone app version
- Language
- Manufacturer
- Phone Model
- Imei

## Update Message

- On update message (currently used for GoLivePhone app updates)

## Activity Monitoring

### Settings

- Daytime hours
- Inactivity alert frequency
- Android ID (for device migration & settings backup)
- Fall detection state (enable/disable)
- Fall detection sensitivity (low/medium/high)
- Fall risk state (enable/disable)
- Inactivity detection state (enable/disable)
- Activity Monitoring state (enable/disable)
- Naps (enable/disable)
- Nap time
- Current phone usage (Belt or pocket)
- Sensor source (Phone or GoLiveClip)

## Logs

- Hourly Records
- Daily Records
- Weekly Records
- Monthly Records
- Inactivity Alerts
- Fall Risk Report

## Location Monitoring

### Settings

- Caretaker alerts (enable/disable)
- Location monitoring state (enable/disable)
- Safe zone monitoring state (enable/disable)
- Safe zone

### Logs

- Location info (provider, latitude, longitude and accuracy)
- Points of interest
- Check ins
- Check outs

### Call Settings

- Android ID (different if more than one device)

### Dialer Settings

- Android ID (different if more than one device)

### Emergency Settings

- Android ID
- Emergency contact (112 or emergency contact)

### Album Settings

- Android ID
- IsSetup
- last sync timestamp
- Media sync (Disabled, 3g or wifi only)

### GoLiveClip

- GoLiveClip state (enable/disable)
- pairedDevice (GoLiveClip mac address)
- pairedDevice firmware version
- LED notifications state (enable/disable)

### Keyboard

- Sound (enable/disable)
- Suggestions (enable/disable)
- Suggestions Language
- Vibration

## Launcher Settings

- Icons (Small/Big)
- Icon order
- Icon selection
- GoLiveLauncher (enable/disable)

## Network Settings

- Mobile data state (enable/disable)
- Wifi state (enable/disable)

## Personal Information

- Address
- Allergies
- Birthdate
- Email
- Gender
- Height
- Medical conditions
- Medications
- Name
- Phone number
- Number of recent falls
- Weight
- Blood type
- Photo

## Privacy Settings

All privacy settings are yes/no fields.

- Send Fall Alert
- Send Inactivity Alert
- Send Location Alert
- Allow Remote Main Menu Configuration
- Send Phone Problems Alert
- Share Activity Level
- Share Location

## Sounds

- Calls tone
- Falls alert tone
- Guide me tone
- Text tone
- Vibration (enable/disable)
- Reminders tone
- Notification
- GoLiveClip tone

## Surveys

- Title
- Description
- Questions
- Visible (y/n)
- Possible answers
- Icon
- Show score (y/n)
- Alert emergency contacts (y/n)

## Survey Results

- Survey ID
- Score
- Status
- Answers

## Vitals

### Blood Pressure

- Body movement (true/false)
- Diastolic pressure
- Irregular pulse (true/false)
- Mean arterial pressure
- Timestamp
- Pulse rate
- Systolic pressure

### Weight

- Weight
- Timestamp

### Unit Settings

- Android ID
- Unit System (metric/imperial)

### Unlocker Settings

- Android ID
- GLP unlocker (enable/disable)

## Alerts

- Fall
- Safezone
- MyStatus
- Battery Phone
- Battery GoLiveClip
- Help Requests
- Connectivity
- Crash
- Exception
- No User Response
- Location Request
- GoLiveClip Link Loss
- GoLiveClip Link Recovery

## System

- Ping status (time)